## ACCENT-200 TOTAL PROTEIN

Nr kat. 7-236 (PL)

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(EN)

## INTENDED USE

TOTAL PROTEIN reagent is intended to determine quantitatively total protein level in plasma and serum. It is intended to monitoring and as an aid in the diagnosis of clinical conditions associated with abnormal protein concentration (hypo-, hyperproteinemia). TOTAL PROTEIN reagent is intended to use on automatic analysers: ACCENT-200 (II GEN) / BS-200, ACCENT-220S / BS-180, ACCENT S120 / BS-230, ACCENT MC240 / BS-240Pro, ACCENT M320 / BS-360E and BS-120.

It is only for *in vitro* diagnostics, for healthcare professional users.

## SUMMARY<sup>1-5</sup>

Determinations of total protein level in serum and plasma are used to diagnosis and monitoring of clinical conditions associated with abnormal protein concentration. Increased total protein level can be observed in multiple myeloma, Waldenstrom macroglobulinemia, dehydration, hyperimmunoglobulinemia caused by polyclonal or monoclonal gammopathies, acute phase response. Decreased values can be observed in malnutrition, malabsorption, liver diseases (cirrhosis, fulminant failure), kidney diseases, protein-losing gastroenteropathies and acute burns.

## METHOD PRINCIPLE<sup>4</sup>

The method is based on the biuret reaction. Protein forms the coloured complex with cupric ions in alkaline medium. The colour intensity is proportional to the protein concentration.

## REAGENTS

### Package

1-REAGENT 3 x 28.5 ml

### Number of tests

ACCENT-200 (II GEN)	250
ACCENT-220S	250
ACCENT S120	410
ACCENT MC240	410
ACCENT M320	410
BS-120	250

## CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

### 1-REAGENT

sodium-potassium tartrate	30 mmol/l
potassium iodide	30 mmol/l
copper sulfate	12 mmol/l
sodium hydroxide	

## REAGENT STABILITY

The reagent, stored at 2-8°C is stable up to expiry date printed on the package. The reagent stored on board of the analyzer at 2-10°C is stable for: 8 weeks (ACCENT-200 (II GEN)), 10 weeks (ACCENT M320) or 11 weeks (ACCENT S120).

## WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients: 1-REAGENT contains sodium hydroxide and copper (II) sulfate pentahydrate.

## Danger

H314 Causes severe skin burns and eye damage.  
H318 Causes serious eye damage.

H411 Toxic to aquatic life with long lasting effects.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing.

Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P273 Avoid release to the environment.

## SPECIMEN<sup>2, 6, 7</sup>

Serum or plasma.

Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Recommended anticoagulants: EDTA, heparine.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

## PROCEDURE

1-REAGENT is ready to use.

Deionized water is recommended as a reagent blank.

For analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, it is recommended to determine the reagent blank during each calibration. Deionized water should be used as reagent blank. When performing calibration, the task type **Calib+Rgt.Blk** should be selected.

### Actions required:

When performing assays in analyzers: ACCENT-200 (II GEN), ACCENT-220S and BS-120, there is a probability of **cross-contamination** affecting the tests results: TOTAL PROTEIN II GEN – FERRUM II GEN, TOTAL PROTEIN II GEN – MICROALBUMIN. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_001ACCENT-200\_CARRYOVER.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analyzers systems: ACCENT- 200 (II GEN), ACCENT MC240, ACCENT M320 and BS-120 the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionized water should be used as a calibrator 0.

For the calibration of automatic analyzers systems: ACCENT-220S, ACCENT S120, the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every: 8 weeks (ACCENT-200 (II GEN)), 10 weeks (ACCENT M320) or 11 weeks (ACCENT S120).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results.

In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

### REFERENCE VALUES 4.8

serum	g/dl	g/l
premature 1 day	3.6 – 6.0	36 – 60
children 1day – 2 weeks	5.3 – 8.3	53 - 83
15 days – 12 months	4.4 – 7.1	44 - 71
1 – 3 years	6.1 – 7.5	61 – 75
3 – 5 years	6.3 – 8.1	63 – 81
6 – 19 years	6.8 – 8.2	68 – 82
adults 20 – 29 years	6.5 – 8.3	65 – 83
30 – 79 years	6.5 – 7.8	65 - 78

Results obtained from plasma analysis might be slightly elevated due to fibrinogen presence:

Sample	Protein concentration increase from serum to plasma
Blood donors	+ 0.25 g/dl
Nonhospitalized patients	+ 0.36 g/dl
Hospitalized patients	+ 0.46 g/dl
Hospitalized patients with CRP<50 mg/dl	+ 0.66 g/dl

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzers ACCENT-200 (II GEN), ACCENT MC240 and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
0.00 g/dl (0.0 g/l) - BS-400
- LoD (Limit of Detection):**  
0.01 g/dl (0.1 g/l) - BS-400
- LoQ (Limit of Quantitation):**  
0.15 g/dl (1.5 g/l) – ACCENT-200 (II GEN)  
0.15 g/dl (1.5 g/l) – ACCENT MC240
- Linearity**  
up to 16.0 g/dl (160 g/l) - ACCENT-200 (II GEN)  
up to 20.0 g/dl (200 g/l) – ACCENT MC240

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Measurement range:**  
0.15 g/dl (1.5 g/l) - 16.0 g/dl (160 g/l) – ACCENT-200 (II GEN)  
0.15 g/dl (1.5 g/l) - 20.0 g/dl (200 g/l) – ACCENT MC240

- Specificity / Interferences**  
Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 950 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [g/dl]	SD [g/dl]	CV [%]
ACCENT-200 (II GEN) n=20	level 1	4.61	0.04	0.82
	level 2	6.85	0.03	0.47
ACCENT MC240 n=20	level 1	4.61	0.02	0.39
	level 2	7.09	0.04	0.56
Reproducibility (day to day)		Mean [g/dl]	SD [g/dl]	CV [%]
ACCENT-200 (II GEN) n=80	level 1	6.9	0.10	1.4
	level 2	4.6	0.04	1.0
ACCENT MC240 n=80	level 1	6.64	0.12	1.8
	level 2	4.37	0.09	2.0

### Method comparison

A comparison between total protein values for samples determined at **ACCENT-200 (II GEN)** (y) and obtained on **BECKMAN COULTER AU680** (x) using 30 plasma samples gave following results:

$$y = 0.9061x + 0.8013 \text{ g/dl};$$
$$R = 0.977 \quad (R - \text{correlation coefficient})$$

A comparison between total protein values for samples determined at **ACCENT-200 (II GEN)** (y) and obtained on **BECKMAN COULTER AU680** (x) using 65 serum samples gave following results:

$$y = 1.0267x - 0.0063 \text{ g/dl};$$
$$R = 0.986 \quad (R - \text{correlation coefficient})$$

A comparison between total protein values for samples determined at **ACCENT MC240** (y) and obtained on **BECKMAN COULTER AU680** (x) using 60 serum samples gave following results:

$$y = 1.0715x - 0.1806 \text{ g/dl};$$
$$R = 0.994 \quad (R - \text{correlation coefficient})$$

A comparison between total protein values for samples determined at **ACCENT MC240** (y) and obtained on **BECKMAN COULTER AU680** (x) using 30 plasma samples gave following results:

$$y = 0.9532x + 0.6738 \text{ g/dl};$$
$$R = 0.986 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT 9

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents:  
1-REAGENT - 18 01 06\*
- Empty packages:  
1-REAGENT - 15 01 10\*
- Wastewater from the analyzer:  
1-REAGENT - 18 01 03\*

### INCIDENTS 10

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

### LITERATURE

- Szczeklik A., Gajewski P., ed. Interna Szczeklika. Mały podręcznik 2018/2019 10th ed., Kraków, Medycyna Praktyczna, (2018)
- Alan H.B., Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B Saunders Company, 916 (2006)
- Chowdhury M., Akhter N., Haque M., Aziz R., Nahar N., Serum Total Protein and Albumin Levels in Different Grades of Protein Energy Malnutrition, *Journal of Bangladesh Society of Physiologist*, 3, 58-60 (2009).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 396, 1780 (2018)
- McPherson R.A., Pincus M.R., ed. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22<sup>nd</sup> ed., ELSEVIER, Philadelphia, (2011).
- WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2 Jan 2002.
- NCCLS - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition, NCCLS document H18-A3 (ISBN 1-56238-555-0), NCCLS, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, USA (2004)
- Thomas L.:Clinical Laboratory Diagnostics. 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 644-647 (1998).

9. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.

10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

### LIST OF CHANGES

Previous version: 04	Current version: 05
Sections updated: <i>INTENDED USE; SUMMARY; CONCENTRATION OF THE ACTIVE INGREDIENTS; WARNING AND NOTES; SPECIMEN; QUALITY CONTROL; REFERENCE VALUES; PERFORMANCE CHARACTERISTICS; WASTE MANAGEMENT; LITERATURE; APPLICATION</i>	
Sections added: <i>REAGENT STABILITY, INCIDENTS</i>	

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- Szczeklik A., Gajewski P., ed. Interna Szczeklika. Mały podręcznik 2018/2019 10th ed., Kraków, Medycyna Praktyczna, (2018)
- Alan H.B., Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B Saunders Company, 916 (2006)
- Chowdhury M., Akhter N., Haque M., Aziz R., Nahar N., Serum Total Protein and Albumin Levels in Different Grades of Protein Energy Malnutrition, *Journal of Bangladesh Society of Physiologist*, 3, 58-60 (2009).
- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 396, 1780 (2018)
- McPherson R.A., Pincus M.R., ed. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22<sup>nd</sup> ed., ELSEVIER, Philadelphia, (2011).
- WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
- NCCLS - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition, NCCLS document H18-A3 (ISBN 1-56238-555-0), NCCLS, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, USA (2004)
- Thomas L.: Clinical Laboratory Diagnostics. 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 644-647 (1998).
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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**Calibration Rule**

Rule	Two-point Linear
Sensitivity	1
Replicates	2
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

**ACCENT- 220S**

**Parameters**

Test	TP	R1	300
No	5	R2	
Full Name	Total Protein	Sample Volume	6
Standard No	5	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	546 nm	Linearity Range	0.2   18.5
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1   19	Factor	
Incuba. Time		<input type="checkbox"/> Prozone check	
Unit	g/dl	q1	<input type="checkbox"/>
Precision	0.01	q2	<input type="checkbox"/>
		q3	<input type="checkbox"/>
		q4	<input type="checkbox"/>
		PC	<input type="checkbox"/>
		Abs	<input type="checkbox"/>

**Calibration Rule**

Rule	Multi-point Linear
Sensitivity	1
Replicates	3
Interval (day)	56
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

**BS-120**

**Parameters**

Test	TP	R1	300
No	5	R2	0
Full Name	Total Protein	Sample Volume	5
Standard No	5	R1 Blank	
Reac. Type	Endpoint	Mixed Rtg. Blank	
Pri. Wave.	546 nm	Linearity Range	0.15   16
Sec. Wave.	670 nm	Linearity Limit	
Direction	Increase	Substrate Limit	
Reac. Time	-1   25	Factor	
Incuba. Time	16	<input type="checkbox"/> Prozone check	
Unit	g/dl	q1	<input type="checkbox"/>
Precision	0.01	q2	<input type="checkbox"/>
		q3	<input type="checkbox"/>
		q4	<input type="checkbox"/>
		PC	<input type="checkbox"/>
		Abs	<input type="checkbox"/>

**Calibration Rule**

Rule	Two-point Linear
Sensitivity	1
Replicates	2
Interval (day)	77
Difference Limit	0
SD	0
Blank Response	0   50000
Error Limit	0
Coefficient	0

● **ACCENT S120**

Chem	TP	No.	005	Sample Type	SERUM
Chemistry	TOTAL PROTEIN	Print name	TOTAL PROTEIN	Reaction Type	Endpoint
Reaction Type	Endpoint	Reaction Direction	positive	Pri Wave	546 nm
Pri Wave	546 nm	Sec Wave	670 nm	Unit	g/dL
Unit	g/dL	Decimal	0.01	Incubation Time	0
Blank Time	-3	-1	Reaction Time	25	27
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol	
Decreased	4	μL	μL	R1	180 μL
Increased	4	μL	20	R2	μL
	μL	μL	180		μL
	μL	μL			μL
	Sample Blank	V	Auto Rerun		
Linearity range (Standard)	0.13	18.00	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-40000	40000
R1 Blank Abs	-40000	40000	On-board Stability	77	Day(s)
Blank Response	-40000	40000	Reagent Alarm Limit		
Twin Chemistry			Enzyme Linear Extension		
Prozone Check					
Q1		Q2	V1	Q3	Q4
Q5		Q6	V3	PC1	PC2
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
		Pretreat Sample Vol		Pretreat Sample Vol	
		μL		μL	
<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>			
Math model	Multi-point linear		Bottle Changed		
Factor			Lot Changed		
Replicates	2		Cal Time		
<b>ACCEPTANCE LIMITS</b>					
Cal Time	1848	Hour			
Slope Diff		SD			
Sensitivity		Repeatability	40000		
Deter Coeff					

● **ACCENT MC240**

Chem	TP	No.	005	Sample Type	SERUM
Chemistry	TOTAL PROTEIN	Print name	TOTAL PROTEIN	Reaction Type	Endpoint
Reaction Type	Endpoint	Reaction Direction	positive	Pri Wave	546 nm
Pri Wave	546 nm	Sec Wave	660 nm	Unit	g/dL
Unit	g/dL	Decimal	0.01	Incubation Time	0
Blank Time	-3	-1	Reaction Time	28	30
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol	
Decreased	4	μL	μL	R1	180 μL
Increased	4	μL	20	R2	μL
	μL	μL	180		μL
	μL	μL			μL
	Sample Blank	V	Auto Rerun		
Linearity range (Standard)	0.15	20.00	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-35000	35000
R1 Blank Abs	-35000	35000	On-board Stability	70	Day(s)
Blank Response	-35000	35000	Reagent Alarm Limit		
Twin Chemistry			Enzyme Linear Extension		
Prozone Check					
Q1		Q2	V1	Q3	Q4
Q5		Q6	V3	PC1	PC2
Sample Pretreatment		Control Pretreatment		Calibrator Pretreatment	
		Pretreat Sample Vol		Pretreat Sample Vol	
		μL		μL	
<b>CALIBRATION SETTINGS</b>		<b>AUTO CALIBRATION</b>			
Math model	Two-point linear		Bottle Changed		
Factor		Replicates	2		
			Lot Changed		
			Cal Time		
<b>ACCEPTANCE LIMITS</b>					
Cal Time	1680	Hour			
Slope Diff		SD			
Sensitivity		Repeatability	35000		
Deter Coeff					

• ACCENT M320

Chem	TP		No.	005		Sample Type	SERUM	
Chemistry	TOTAL PROTEIN		Print name	TOTAL PROTEIN		Reaction Direction	positive	
Reaction Type	Endpoint		Sec Wave	660 nm		Decimal	0.01	
Pri Wave	546 nm		Incubation Time	0		Reaction Time	28 30	
Unit	g/dL		Sample Vol	4		Aspirated	20	
Blank Time	-3 -1		Diluent	180		Reagent Vol	180	
Standard	4		Decreased	4		R1	180	
Increased			Sample Blank	<input type="checkbox"/>		Auto Rerun	<input checked="" type="checkbox"/>	

Linearity range (Standard)	0.13 18.00		Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-35000 35000	
R1 Blank Abs	-35000 35000		On-board Stability	70 Dav(s)	
Blank Response	-35000 35000		Reagent Alarm Limit		
Twin Chemistry	<input type="checkbox"/> Prozone Check		<input type="checkbox"/> Enzyme Linear Extension		
Q1	Q2	V1	O3	Q4	V2
Q5	Q6	V3	PC1	PC2	
<input type="checkbox"/> Sample Pretreatment	<input type="checkbox"/> Control Pretreatment	<input type="checkbox"/> Calibrator Pretreatment	Pretreat Sample Vol	μL	

**CALIBRATION SETTINGS**

Math model	Two-point linear		<input type="checkbox"/> Bottle Changed
Factor	Replicates	2	

**AUTO CALIBRATION**

<input type="checkbox"/> Lot Changed	<input type="checkbox"/> Cal Time
--------------------------------------	-----------------------------------

**ACCEPTANCE LIMITS**

Cal Time	1680 Hour	
Slope Diff	SD	
Sensitivity	Repeatability	35000
Deter Coeff		

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