

# CORMAY IMMUNO-MULTICAL

KARTA WARTOŚCI WYZNACZONYCH / ASSAY SHEET



F07-QAD-01\_06 15.12.2023

<b>Nr kat. / Cat. No</b>	<b>Seria / Lot</b>	<b>Data ważności / Expiration date</b>
<b>4-287</b>	<b>312-72</b>	<b>2025-07</b>

PARAMETR / PARAMETER	IMMUNO-MULTICAL				
	1	2	3	4	5
<b><math>\alpha</math>-1 ANTYTRYPSYNA / <math>\alpha</math>-1 ANTITRYPSIN</b>	0.17 g/L	0.34 g/L	0.68 g/L	1.39 g/L	2.72 g/L
<b>KWAŚNA <math>\alpha</math>1-GLIKOPROTEINA / <math>\alpha</math>-1 GLYCOPROTEIN ACID</b>	0.13 g/L	0.25 g/L	0.50 g/L	1.00 g/L	2.00 g/L
<b>ANTITHROMBIN III/ ANTYTROMBINA III</b>	0.038 g/L	0.076 g/L	0.147 g/L	0.302 g/L	0.609 g/L
<b>SKŁADNIK C3 DOPELNIACZA / COMPLEMENT C3</b>	0.23 g/L	0.47 g/L	0.92 g/L	1.87 g/L	3.77 g/L
<b>SKŁADNIK C4 DOPELNIACZA / COMPLEMENT C4</b>	0.044 g/L	0.088 g/L	0.175 g/L	0.350 g/L	0.705 g/L
<b>CERULOPLAZMINA / CERULOPLASMIN</b>	0.082 g/L	0.142 g/L	0.264 g/L	0.533 g/L	1.011 g/L
<b>HAPTOGLOBINA / HAPTOGLOBIN</b>	0.17 g/L	0.35 g/L	0.70 g/L	1.41 g/L	2.84 g/L
<b>IgA</b>	0.32 g/L	0.67 g/L	1.33 g/L	2.67 g/L	5.37 g/L
<b>IgG</b>	1.62 g/L	3.26 g/L	6.51 g/L	13.06 g/L	26.28 g/L
<b>IgM</b>	0.14 g/L	0.27 g/L	0.55 g/L	1.10 g/L	2.21 g/L
<b>TRANSFERYNA /TRANSFERRIN</b>	0.42 g/L	0.84 g/L	1.66 g/L	3.38 g/L	6.81 g/L