

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 304-74

Data ważn. / Exp.: 2025-01

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	0.89 mg/dl	0.62 – 1.16 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	146 IU/ml	117 – 175 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	16.9 IU/ml	12.68 – 21.13 IU/ml

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 304-74

Data ważn. / Exp.: 2025.01

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	0.89 mg/dl	0.62 – 1.16 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	146 IU/ml	117 – 175 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	16.9 IU/ml	12.68 – 21.13 IU/ml

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 304-72

Data ważn. / Exp.: 2025.01

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	2.44 mg/dl	2.07 – 2.81 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	293 IU/ml	249 – 337 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	35.5 IU/ml	26.63 – 44.38 IU/ml

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 304-72

Data ważn. / Exp.: 2025.01

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	2.44 mg/dl	2.07 – 2.81 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	293 IU/ml	249 – 337 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	35.5 IU/ml	26.63 – 44.38 IU/ml