

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 307-76

Data ważn. / Exp.: 2025-05

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	0.97 mg/dl	0.68 – 1.26 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	150 IU/ml	120 – 180 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	16.1 IU/ml	12.08 – 20.13 IU/mL

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 307-76

Data ważn. / Exp.: 2025-05

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	0.97 mg/dl	0.68 – 1.26 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	150 IU/ml	120 – 180 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	16.1 IU/ml	12.08 – 20.13 IU/mL

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 307-74

Data ważn. / Exp.: 2025-05

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	2.57 mg/dl	2.18 - 2.96 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	310 IU/ml	264 – 357 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	37.1 IU/ml	27.83 – 46.38 IU/ml

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 307-74

Data ważn. / Exp.: 2025-05

PARAMETR / COMPONENT	METODA / METHOD	WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE	ZAKRES / RANGE
BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP)	TIA-wzmacniane lateksem Latex-enhanced TIA	2.57 mg/dl	2.18 - 2.96 mg/dl
ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO)	TIA-wzmacniane lateksem Latex-enhanced TIA	310 IU/ml	264 – 357 IU/ml
CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF)	TIA-wzmacniane lateksem Latex-enhanced TIA	37.1 IU/ml	27.83 – 46.38 IU/ml