

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 212-73

Data ważn. / Exp.: 2024.10

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|-------------------|
| BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP) | TIA-wzmacniane lateksem Latex-enhanced TIA | 0.90 mg/dl | 0.63 – 1.17 mg/dl |
| ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO) | TIA-wzmacniane lateksem Latex-enhanced TIA | 141 IU/ml | 113 – 169 IU/ml |
| CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF) | TIA-wzmacniane lateksem Latex-enhanced TIA | 18.4 IU/ml | 13.8 – 23.0 IU/ml |

CORMAY IMMUNO-CONTROL I (L)

Nr kat. / Cat. No 4-288

Seria / Lot: 212-73

Data ważn. / Exp.: 2024.10

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|-------------------|
| BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP) | TIA-wzmacniane lateksem Latex-enhanced TIA | 0.90 mg/dl | 0.63 – 1.17 mg/dl |
| ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO) | TIA-wzmacniane lateksem Latex-enhanced TIA | 141 IU/ml | 113 – 169 IU/ml |
| CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF) | TIA-wzmacniane lateksem Latex-enhanced TIA | 18.4 IU/ml | 13.8 – 23.0 IU/ml |

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 212-71

Data ważn. / Exp.: 2024.10

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|---------------------|
| BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP) | TIA-wzmacniane lateksem Latex-enhanced TIA | 2.52 mg/dl | 2.14 – 2.90 mg/dl |
| ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO) | TIA-wzmacniane lateksem Latex-enhanced TIA | 293 IU/ml | 249 – 337 IU/ml |
| CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF) | TIA-wzmacniane lateksem Latex-enhanced TIA | 37.3 IU/ml | 27.98 – 46.63 IU/ml |

CORMAY IMMUNO-CONTROL I (H)

Nr kat. / Cat. No 4-288

Seria / Lot: 212-71

Data ważn. / Exp.: 2024.10

| PARAMETR / COMPONENT | METODA / METHOD | WARTOŚĆ ŚRODKOWA / ASSIGNED VALUE | ZAKRES / RANGE |
|---|---|-----------------------------------|---------------------|
| BIALKO C-REAKTYWNE (CRP) C-REACTIVE PROTEIN (CRP) | TIA-wzmacniane lateksem Latex-enhanced TIA | 2.52 mg/dl | 2.14 – 2.90 mg/dl |
| ANTY-STREPTOLIZYNA O (ASO) ANTY-STREPTOLYSIN O (ASO) | TIA-wzmacniane lateksem Latex-enhanced TIA | 293 IU/ml | 249 – 337 IU/ml |
| CZYNNIK REUMATOIDALNY (RF) RHEUMATOID FACTOR (RF) | TIA-wzmacniane lateksem Latex-enhanced TIA | 37.3 IU/ml | 27.98 – 46.63 IU/ml |