




G-CAL

CALIBRATOR

LOT GC0326

 2026-06-30

Parameters	Mythic 22		Mythic 22 AL		Mythic 60		Mythic 70	
	<i>Mean</i>	<i>Range</i>	<i>Mean</i>	<i>Range</i>	<i>Mean</i>	<i>Range</i>	<i>Mean</i>	<i>Range</i>
FSC Target	84	83-85	82	81-83	79	78-80	78	77-79
ALL Target	63	62-64	63	62-64	66	65-67	66	65-67

ORPHEE SA

19, chemin du champ-des-filles
CH-1228 Geneva / Plan-les-Ouates
SWITZERLAND