

**HAEM 8 CONTROL L**  
**HAEM 8 CONTROL N**  
**HAEM 8 CONTROL H**



KARTA WARTOŚCI WYZNACZONYCH / ASSAY SHEET / АНАЛИТИЧЕСКА ЛИСТА

F07-QAD-01\_15 15.12.2023

Nr kat. / Cat. No / Номер кат.	Seria / Lot / Серия	Data ważności / Expiration date / Срок годности
<b>8-601</b>	<b>403-X0424</b>	<b>2024.07.05</b>
<b>8-602</b>		
<b>8-603</b>		

**Mindray BC 3000 PLUS**

(odczynniki / reagents / реагенты PZ CORMAY S.A.)

Parametr / Parameter / Параметр		HAEM 8 CONTROL					
		LOW		NORMAL		HIGH	
		Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±
WBC	10 <sup>9</sup> /l	2.1	± 0.5	10.0	± 1.0	20.0	± 2.5
RBC	10 <sup>12</sup> /l	2.32	± 0.18	4.71	± 0.24	5.70	± 0.30
Hb	g/dl	6.3	± 0.4	13.5	± 0.6	18.1	± 0.8
Hct	%	17.6	± 1.5	40.1	± 2.0	54.6	± 2.4
MCV	fl	75.9	± 5.0	85.1	± 5.0	95.8	± 5.0
MCH	pg	27.2	± 2.5	28.8	± 2.5	31.8	± 2.5
MCHC	g/dl	35.8	± 3.0	33.8	± 3.0	33.2	± 3.0
Plts	10 <sup>9</sup> /l	66	± 20	239	± 40	490	± 60
MPV	fl	8.4	± 3.0	8.2	± 3.0	8.2	± 3.0

**Mythic 18 Orphee**

(odczynniki / reagents / реагенты PZ CORMAY S.A.)

Parametr / Parameter / Параметр		HAEM 8 CONTROL					
		LOW		NORMAL		HIGH	
		Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±
WBC	10 <sup>9</sup> /l	2.3	± 0.5	10.3	± 1.0	20.9	± 2.5
RBC	10 <sup>12</sup> /l	2.25	± 0.18	4.51	± 0.24	5.47	± 0.30
Hb	g/dl	6.2	± 0.4	13.0	± 0.6	17.3	± 0.8
Hct	%	17.2	± 1.5	37.8	± 2.0	51.5	± 2.4
MCV	fl	76.6	± 5.0	83.9	± 5.0	94.1	± 5.0
MCH	pg	27.5	± 2.5	28.8	± 2.5	31.7	± 2.5
MCHC	g/dl	35.9	± 3.0	34.4	± 3.0	33.7	± 3.0
Plts	10 <sup>9</sup> /l	74	± 20	271	± 40	536	± 60
MPV	fl	8.5	± 3.0	8.1	± 3.0	8.2	± 3.0

**ABX: Micros (CBC-mode)**

(odczynniki / reagents / реагенты PZ CORMAY S.A.)

Parametr / Parameter / Параметр		HAEM 8 CONTROL					
		LOW		NORMAL		HIGH	
		Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±	Średnia / Mean / Средняя	Zakres ± / Range ± / Предел ±
WBC	10 <sup>9</sup> /l	2.3	± 0.5	10.6	± 1.0	22.0	± 2.5
RBC	10 <sup>12</sup> /l	2.21	± 0.18	4.55	± 0.24	5.55	± 0.30
Hb	g/dl	6.2	± 0.4	13.3	± 0.6	17.5	± 0.8
Hct	%	16.0	± 1.5	37.2	± 2.0	51.4	± 2.4
MCV	fl	72.0	± 5.0	82.0	± 5.0	93.0	± 5.0
MCH	pg	28.3	± 2.5	29.1	± 2.5	31.5	± 2.5
MCHC	g/dl	39.1	± 3.0	35.6	± 3.0	34.0	± 3.0
Plts	10 <sup>9</sup> /l	69	± 20	253	± 40	494	± 60
MPV	fl	9.4	± 3.0	8.8	± 3.0	8.7	± 3.0