



# MYT-3CRP

## HEMATOLOGY CONTROLS

### CONTROL

**LOT** B0224



2024-05-10

### ASSAY VALUES AND EXPECTED RANGES

Instrument	Parameter	Low		Normal		High	
		<b>LOT</b> B0224L	B0224L	<b>LOT</b> B0224N	B0224N	<b>LOT</b> B0224H	B0224H
Orphée Mythic 3CRP	WBC × 10 <sup>3</sup> /μL		2.04 ± 0.60		8.19 ± 1.00		21.73 ± 2.50
	LYM %		67.9 ± 9.0		35.3 ± 8.0		14.5 ± 8.0
	GRAN %		29.5 ± 9.0		59.6 ± 8.0		79.3 ± 8.0
	MID %		2.6 ± 2.6		5.1 ± 5.1		6.2 ± 4.0
	LYM × 10 <sup>3</sup> /μL		1.39 ± 0.30		2.89 ± 0.70		3.15 ± 2.00
	GRAN × 10 <sup>3</sup> /μL		0.60 ± 0.30		4.88 ± 0.70		17.23 ± 2.00
	MID × 10 <sup>3</sup> /μL		0.05 ± 0.05		0.42 ± 0.42		1.35 ± 0.50
	RBC × 10 <sup>6</sup> /μL		2.41 ± 0.24		4.65 ± 0.24		5.71 ± 0.30
	HGB g/dL		58 ± 6		136 ± 6		187 ± 8
	HCT %		18.7 ± 2.0		41.3 ± 3.0		54.7 ± 4.0
	MCV fL		77.5 ± 5.0		88.8 ± 5.0		95.8 ± 5.0
	MCH pg		24.0 ± 2.5		29.0 ± 2.5		32.5 ± 2.5
	MCHC g/dL		317 ± 30		331 ± 30		347 ± 30
	RDW-CV %		18.0 ± 3.0		14.7 ± 3.0		14.3 ± 3.0
	RDW-SD fL		51.2 ± 8.0		48.6 ± 8.0		50.5 ± 8.0
	PLT × 10 <sup>3</sup> /μL		71 ± 20		272 ± 40		549 ± 60
	MPV fL		8.9 ± 3.0		9.4 ± 3.0		9.4 ± 3.0
	PDW fL		11.6 ± 3.0		12.4 ± 3.0		12.7 ± 3.0
	PCT %		0.063 ± 0.063		0.254 ± 0.200		0.512 ± 0.200
	P-LCR %		24.2 ± 10.0		28.1 ± 10.0		28.2 ± 10.0
P-LCC × 10 <sup>3</sup> /μL		17 ± 15		76 ± 25		153 ± 30	

**NOTE** Before using, refer to the instruction sheet for mixing directions.  
All brands and products are trademarks or registered trademarks of their respective companies.



B0224L



B0224N



B0224H

**ORPHÉE SA**  
19 chemin du Champ-des-Filles  
CH-1228 Plan-les-Ouates / Geneva  
Switzerland

**PZ CORMAY S.A.**  
ul. Wiosenna 22,  
05-092 Łomianki,  
POLSKA

