



# MYT-3CRP

## HEMATOLOGY CONTROLS

### CONTROL

**LOT** B0225

2025-05-10

#### ASSAY VALUES AND EXPECTED RANGES

Instrument	Parameter	Low		Normal		High	
		<b>LOT</b>	B0225L	<b>LOT</b>	B0225N	<b>LOT</b>	B0225H
Orphée Mythic 3CRP	WBC × 10 <sup>3</sup> /μL		2.00 ± 0.60		7.90 ± 1.00		21.30 ± 2.50
	LYM %		67.6 ± 9.0		35.3 ± 8.0		12.9 ± 8.0
	GRAN %		29.6 ± 9.0		62.0 ± 8.0		84.1 ± 8.0
	MID %		2.8 ± 2.8		2.7 ± 2.7		3.0 ± 3.0
	LYM × 10 <sup>3</sup> /μL		1.35 ± 0.30		2.79 ± 0.70		2.75 ± 2.00
	GRAN × 10 <sup>3</sup> /μL		0.59 ± 0.30		4.90 ± 0.70		17.91 ± 2.00
	MID × 10 <sup>3</sup> /μL		0.06 ± 0.06		0.21 ± 0.21		0.64 ± 0.64
	RBC × 10 <sup>6</sup> /μL		2.43 ± 0.24		4.71 ± 0.24		5.56 ± 0.30
	HGB g/dL		5.8 ± 0.6		13.7 ± 0.6		18.2 ± 0.8
	HCT %		19.2 ± 2.0		42.1 ± 3.0		54.0 ± 4.0
	MCV fL		79.1 ± 6.0		89.4 ± 6.0		97.2 ± 7.0
	MCH pg		24.1 ± 2.5		29.4 ± 2.5		33.0 ± 2.5
	MCHC g/dL		31.1 ± 3.0		33.3 ± 3.0		34.6 ± 3.0
	RDW-CV %		17.2 ± 3.0		15.1 ± 3.0		14.4 ± 3.0
	RDW-SD fL		50.2 ± 8.0		50.1 ± 8.0		51.3 ± 8.0
	PLT × 10 <sup>3</sup> /μL		75 ± 20		260 ± 40		537 ± 60
	MPV fL		9.0 ± 3.0		9.2 ± 3.0		8.9 ± 3.0
	PDW fL		11.7 ± 3.0		11.8 ± 3.0		11.3 ± 3.0
	PCT %		0.066 ± 0.066		0.231 ± 0.200		0.461 ± 0.200
	P-LCR %		25.1 ± 10.0		25.8 ± 10.0		23.7 ± 10.0
P-LCC × 10 <sup>3</sup> /μL		18 ± 15		65 ± 25		124 ± 30	

**NOTE** Before using, refer to the instruction sheet for mixing directions.  
All brands and products are trademarks or registered trademarks of their respective companies.



B0225L



B0225N



B0225H

**ORPHÉE SA**  
19 chemin du Champ-des-Filles  
CH-1228 Plan-les-Ouates / Geneva  
Switzerland



**PZ CORMAY S.A.**  
ul. Wiosenna 22,  
05-092 Łomianki,  
POLSKA

