



MYT-3CRP

HEMATOLOGY CONTROLS

CONTROL

LOT B1124



2025-02-10

ASSAY VALUES AND EXPECTED RANGES

Instrument	Parameter	Low		Normal		High	
		LOT B1124L	B1124L	LOT B1124N	B1124N	LOT B1124H	B1124H
Orphée Mythic 3CRP	WBC × 10 ³ /μL		1.96 ± 0.60		7.76 ± 1.00		20.67 ± 2.50
	LYM %		68.3 ± 9.0		35.6 ± 8.0		17.6 ± 8.0
	GRAN %		28.7 ± 9.0		59.4 ± 8.0		76.9 ± 8.0
	MID %		3.0 ± 3.0		5.0 ± 5.0		5.5 ± 5.0
	LYM × 10 ³ /μL		1.34 ± 0.30		2.76 ± 0.70		3.63 ± 2.00
	GRAN × 10 ³ /μL		0.56 ± 0.30		4.61 ± 0.70		15.90 ± 2.00
	MID × 10 ³ /μL		0.06 ± 0.06		0.39 ± 0.39		1.14 ± 1.00
	RBC × 10 ⁶ /μL		2.40 ± 0.24		4.63 ± 0.24		5.64 ± 0.30
	HGB g/dL		5.8 ± 0.6		13.4 ± 0.6		18.6 ± 0.8
	HCT %		18.8 ± 2.0		40.5 ± 3.0		54.4 ± 4.0
	MCV fL		78.4 ± 5.0		87.5 ± 5.0		96.5 ± 5.0
	MCH pg		24.5 ± 2.5		29.5 ± 2.5		33.6 ± 2.5
	MCHC g/dL		31.7 ± 3.0		33.9 ± 3.0		35.3 ± 3.0
	RDW-CV %		16.7 ± 3.0		15.6 ± 3.0		14.3 ± 3.0
	RDW-SD fL		48.7 ± 8.0		51.1 ± 8.0		51.1 ± 8.0
	PLT × 10 ³ /μL		67 ± 20		260 ± 40		530 ± 60
	MPV fL		9.1 ± 3.0		9.4 ± 3.0		9.2 ± 3.0
	PDW fL		11.7 ± 3.0		12.2 ± 3.0		12.2 ± 3.0
	PCT %		0.063 ± 0.063		0.250 ± 0.200		0.498 ± 0.200
	P-LCR %		25.7 ± 10.0		27.5 ± 10.0		26.6 ± 10.0
P-LCC × 10 ³ /μL		18 ± 15		73 ± 25		144 ± 30	

NOTE Before using, refer to the instruction sheet for mixing directions.
All brands and products are trademarks or registered trademarks of their respective companies.



B1124L



B1124N



B1124H

ORPHÉE SA
19 chemin du Champ-des-Filles
CH-1228 Plan-les-Ouates / Geneva
Switzerland



PZ CORMAY S.A.
ul. Wiosenna 22,
05-092 Łomianki,
POLSKA

