



# MYTHIC CRP Controls

CRP CONTROLS

**CONTROL**

**LOT** 2023121301

 2024-12-12

## ASSAY VALUES AND EXPECTED RANGES

Instrument	Control	Reference Value (mg/L)	Deviation Range (mg/L)
Orphée Mythic 3CRP	Level 1	3.3	± 1.0
	Level 2	14.8	± 4.4
	Level 3	41.7	± 12.5

**NOTE** Before using, refer to the instruction sheet for mixing directions.

All brands and products are trademarks or registered trademarks of their respective companies.

### ORPHÉE SA

19 chemin du Champ-des-Filles  
CH-1228 Plan-les-Ouates / Geneva  
Switzerland



### PZ CORMAY S.A.

ul. Wiosenna 22,  
05-092 Łomianki, POLSKA

