

Clinical Chemistry

Intelligent Technology
You Can Trust




The lifeline to better health

We are constantly developing innovative In-Vitro Diagnostic solutions to accurately uncover the truth hidden in a single drop of blood. Our actions help to understand the known and discover the unknown, leading to the earliest detection, treatment and monitoring of the health of patients all over the world.

Our cutting-edge In-Vitro Diagnostic solutions provide the lifeline that connects doctors and clinicians to their patients to deliver a positive health outcome.

For Almost 40 Years We Connect Healthcare Professionals with European Quality & Innovation in IVD

Cormay Diagnostics has a lot of experience and know-how in producing comprehensive range of market leading clinical chemistry and hematology reagents, together with innovations and patented technology have allowed us to deliver outstanding performance and reliability of our IVD solutions.
Our Cormay and Orphée IVD brands have become trusted by healthcare professionals from all over the world.



Global Range

We act in over 120 countries and connect distributors all over the world with world class technology & performance for small & medium size clinical laboratories.

We act in
120
countries
worldwide

Over
20k
analyzers
installed
worldwide

Over
350k
kits produced
per year



We connect people in need with tests that
can help them in their fight for health.

Our unwavering commitment to providing healthcare professionals with more effective tests that guarantee consistent, high - quality results and expert diagnosis ensures that patients around the world can count on a quick and healthy recovery.



The Accent Family

Of Intelligent Multi-Patented Clinical Chemistry Analyzers

Providing the **BEST POSSIBLE PATIENT CARE** is our top priority. We understand physicians rely on fast and accurate diagnostic information to prescribe the appropriate course of action for a speedy, healthy recovery. The sooner the problem is defined, the sooner patients can start their treatment with a higher probability of a better outcome, where early diagnosis is crucial to treatment's success.

First launched in 2007, Cormay's Accent clinical chemistry analysers have become the new standard for clinicians having rapidly established a reputation for cost effective, high quality, reliable performance in clinical laboratories all over the world.



UNCOMPROMISED
PATIENT
DIAGNOSIS



UNINTERRUPTED
RUN



ADVANCED
SOFTWARE




SIGNIFICANT
COST SAVINGS



New Generation of fully automated, bench - top Accent analyzers offers technological advancements and intelligence that ensure you work smarter, faster and with greater accuracy.

ACCENT S120

- ✓ Up to 200 T/H* & up to 400 T/H* with ISE
- ✓ 40 Disposable Cuvettes
- ✓ 8 Wavelengths



ACCENT MC240

- ✓ 240 T/H Constant & 400 T/H with ISE
- ✓ 80 Reusable Cuvettes
- ✓ Clog Detection
- ✓ Grating System with 12 Wavelengths



ACCENT M320

- ✓ Up to 360 T/H* & up to 540 T/H* with ISE
- ✓ 80 Reusable Cuvettes
- ✓ Clog Detection
- ✓ Grating System with 12 Wavelengths

* theoretical value

Uncover the Secrets To Your Health

When a diagnosis is accurate and made in a timely manner, a patient has the best opportunity for a positive health outcome because clinical decision making will be tailored to a correct understanding of the patient's health problem.



Market leading
Clinical Chemistry reagents from Europe

48
months
the longest
shelf life

12
weeks
the longest
on board stability

12
weeks
the longest
calibrations stability

Allergic	Anemia	Bone metabolism	Cardiac	Diabetic	Toxicological
Total IgE	Ferritin Iron Transferrin UIBC	ALP Calcium Phosphorus	ALAT Apo B ASAT Cholesterol HDL Direct LDL Direct CK CK-MB CRP Glucose Lactate LDH Lipoprotein (a) Myoglobin Potassium ISE Sodium ISE Triglycerides	Cholesterol HDL Direct LDL Direct Creatinine Cystatin C Glucose HbA1c Direct Lactate Microalbumin Triglycerides	Ethanol

Comprehensive Range of Panels
and Parameters for the Conclusive Diagnosis

Electrolytes	Hepatic	Infectious	Covid - 19	Routine Screening	Inflammatory	Immune	Lipids	Nutritional	Oncological	Pancreatic	Pediatric Screening	Pregnancy	Renal	Specific Proteins
Calcium Chloride ISE Iron Magnesium Phosphorus Potassium ISE Sodium ISE	AAT ALAT Albumin ALP ASAT Bilirubin Direct Bilirubin Total Ceruloplasmin Cholinesterase GGT Lactate Total Bile Acids Total Protein Transferrin	AAT AGP ASO Ceruloplasmin C3 C4 CRP Haptoglobin	ALAT Albumin ASAT Bilirubin Total Creatinine CRP Ferritin LDH Procalcitonin Urea	ALAT Albumin ALP Amylase ASAT Bilirubin Total Bilirubin Direct Calcium Chloride ISE Cholesterol HDL Direct LDL Direct CK Creatinine GGT Glucose Iron LDH Lipase Magnesium Phosphorus Potassium ISE Sodium ISE Triglycerides Total Protein Urea Uric Acid Urine Proteins	AAT AGP ASO Ceruloplasmin C3 C4 CRP Ferritin Haptoglobin Procalcitonin RF	C3 C4 IgA IgG IgM Total IgE	Apo B Cholesterol HDL Direct LDL Direct Lipoprotein (a) Triglycerides	Albumin Iron Magnesium Total Protein Transferrin Uric Acid Urea	AGP Calcium Cholinesterase Ferritin IgA IgG IgM Lactate LDH Lipase Total Bile Acids Total IgE Uric Acid	Amylase Glucose Lipase	Albumin AAT ASO Bilirubin Direct Bilirubin Total CRP Creatinine Cystatin C Glucose IgA IgG IgM Iron Lactate Magnesium Phosphorus Total IgE Total Protein Transferrin	ALAT ALP ASAT Bilirubin Total Calcium Creatinine CRP Ferritin Glucose GGT Iron Magnesium Potassium ISE Sodium ISE Total Bile Acids Transferrin UIBC	Albumin Calcium Ceruloplasmin Creatinine Cystatin C Glucose Lactate Microalbumin Phosphorus Potassium ISE Sodium ISE Urea Uric Acid Urine Proteins	AAT AGP Apo B ASO Ceruloplasmin C3 C4 CRP Cystatin C Ferritin Haptoglobin HbA1c Direct IgA IgG IgM Lipoprotein (a) Microalbumin Myoglobin Procalcitonin RF Total IgE Transferrin

Confidence in our Reagents


- High Quality Raw Material
- Tested and Proven Formulations
- State-of-the-Art Production Technology
- Extensive Quality Control

Diagnostic Test Results You Can Trust

Accent Multi - Patented Performance Enhancement Technology

At Cormay, we are driven by both clinicians and patients fundamental need for fast and accurate results. Our relentless efforts and meticulous attention to detail has resulted in the development of state - of - the - art clinical chemistry analyzer that sets new standards in diagnostic performance.

Accent combines powerful Patented Optical Technology and Smart On - Board Sample Management with Cormay's comprehensive portfolio of market leading clinical chemistry reagents to deliver a diagnostic test results you can trust.



Risk from Carryover Significantly Reduced

Accent's Smart On - Board Management incorporates multiple carryover reducing measures including waterfall technology for internal and external probe cleaning, a washing station optimized for thorough cuvette washing and independent mixing bar.


Carryover reduced to $\leq 0.05\%$ (for sample) and ≤ 1 ppm (for washing station)


Outstanding Accuracy & Linearity

Accent's Patented Concave Grating Light-Splitting Technology sets new standards in clinical chemistry performance.

High precision, high quality optical components deliver enhanced sensitivity, the highest signal - to - noise performance and lowest stray light for a diagnostic test result you can trust.

Enhanced Signal - to - Noise PDA Performance through Simultaneous Background Correction





The Complete Laboratory Solution Intelligent Technology You Can Trust

With today's clinical laboratories increasingly under time and cost pressure to deliver fast and accurate results, Accent provide productivity to small and medium sized laboratories without compromising patient care. Accent's extended test menu, built - in intelligence, smart software and lower reagent consumption brings the performance of top end analysers to your laboratory.

Only with our reagents, calibrators and controls we can make full use of Accent's technology and bring a **powerfull SYSTEM**.

And together with the deep involvement and support of our customers business, we provide the **COMPLETE VALIDATED LABORATORY SOLUTION**.

PRODUCTIVITY



Uninterrupted Run

Critical samples and reagents loading during run without workflow disruption.



Reduced Instrument Downtime

Maximize time for critical testing of patient samples through Accent's Self - Diagnostic Service and Maintenance Software.



Significantly Reduced Reruns

Avoid unnecessary and time wasting work - flow disruptions with Accent's Intelligent Probe Technology.

- ✓ Clog Detection
- ✓ Collision Protection
- ✓ Liquid Level Detection
- ✓ Inventory Monitoring
- ✓ Reagent Pre-Warming



Conclusive Diagnostic Test Results the First Time

Immediately confirm flagged patient results during unattended sample runs with Accent's Intelligent REFLEX option.

Automated detection of interferents (hemolysis, icterus and lipemia) with SERUM INDEX Program helps reduce pre - analytical errors affecting test results.




Cost Savings You Can Count On

Optimized applications which significantly reduce reagent consumption without compromising performance. 100 µl minimum reaction volumes.

Ready to use reagents for a broad number of clinical chemistry analyzers

As a well-known provider of clinical chemistry solutions, Cormay develops and manufactures high quality reagents for open and dedicated lines.

Our reagents fit the broad number analyzers and ENHANCE EACH AND EVERY OF THEM.



User & Environmental Friendly Reagents

During the years, Cormay has started changes in the reagent formulations to find alternative detergents which are more eco-friendly for environment.

From many reagents like: ALAT, ASAT, BILIRUBIN DIRECT, CALCIUM ARSENAZO, LIPASE, PHOSPHORUS, URIC ACID, UREA, we removed Octylphenol Ethoxylates (OPEs), such as Triton X, which have toxic effect on aquatic species and replaced them with environment friendly detergents.

In case of other reagents like Concentrated Detergent B – we made biodegradable detergent formulation and removed phosphate – component which can cause nutrient pollution of water and overgrowth of algae – which can have serious impacts on the environment and human health.

WE constantly change for You and for ENVIRONMENTAL protection!

Join to thousands of
satisfied customers



Just scan the QR code
for more marketing materials




CLINICAL CHEMISTRY



HEMATOLOGY



VETERINARY



IMMUNOLOGY



URINALYSIS



RAPID TESTS



CORMAY



303 Pulawska St.
02-785 Warszawa | Poland
tel.: +48 (22) 751 79 10
international@cormaydiagnostics.com
www.cormaydiagnostics.com