

## ADVISORY NOTE

Date of issue	17.10.2023
Note number	ND_2310_02
Product name and catalogue number	Myt-5D (Cat. No MYT502, MYT506N)
Lot number	MYT2310
Expiry date	2023-12-05
Issue description	
<p>APPLIES TO USERS OF <b>MYTHIC 70</b> ANALYZER!!!</p> <p>During re-verification of the value sheet for Myt-5D control blood it was found that <b>the values for PLT, IMM, ALY and HCT were incorrectly provided.</b></p> <p>To solve the problem, the values were redetermined and reported in revised assay value sheet (Attachment 1).</p>	
Necessary actions	
<p><b>Read carefully the attachment and use the re-determined values of PLT, IMM, ALY and HCT parameters for control blood.</b></p> <p>Please complete and return the verification form.</p> <p>The note is only for reference. <b>Do not discard the product!</b></p>	

Please provide this information to all Users of the above product and keep the Advisory Note as a quality system documentation.

Please complete and return the attached Verification Form to confirm receipt of the Advisory Note.

We apologize for any inconvenience caused by the above situation.

Sincerely,

Quality Assurance Manager