

**ADVISORY NOTE**

Date of issue	29.02.2024
Note number	ND_2402_03
Product name and catalogue numer	Myt-5D (Cat. No MYT502, MYT506N)
Lot number	MYT2402
Expiry date	2024-04-05
Issue description	
<p>APPLIES TO USERS OF <b>MYTHIC 70</b> ANALYZER!!!</p> <p>During re-verification of the value sheet for Myt-5D control blood (level Normal and High) it was found that <b>the values for PLT, MPV, PCT, PDW and PLCC were incorrectly provided.</b></p> <p>To solve the problem, the values were redetermined and reported in revised assay value sheet (Attachment 1).</p>	
Necessary actions	
<p><b>Read carefully the attachment and use the re-determined values of PLT, MPV, PCT, PDW and PLCC parameters for control blood.</b></p> <p>Please complete and return the verification form.</p> <p>The note is only for reference. <b>Do not discard the product!</b></p>	

Please provide this information to all Users of the above product and keep the Advisory Note as a quality system documentation.

Please complete and return the attached Verification Form to confirm receipt of the Advisory Note.

We apologize for any inconvenience caused by the above situation.

Sincerely,

Quality Assurance Manager