







CHANGE NOTIFICATION

Notification number	PZ_2309_02	Date	20.09.2023
Change number	098		
Subject of the change	Classification and labelling change.		
Affected products	Cat no	Product name	
	6-251	CORMAY ASO Latex 100	
	6-254	CORMAY CRP Latex 100	
	6-252	CORMAY RF Latex 100	
	6-256	CORMAY RPR Carbon 500	
	6-259	CORMAY TPHA 100	
Change category	<input type="checkbox"/> raw materials <input type="checkbox"/> quality control <input type="checkbox"/> packages <input checked="" type="checkbox"/> informations provided with the product <input type="checkbox"/> manufacturing <input type="checkbox"/> final product <input type="checkbox"/> other		
Change description	<p>The classification and labelling of agglutination tests has been changed:</p> <ul style="list-style-type: none"> • CORMAY CRP Latex 100, • CORMAY ASO Latex 100 • CORMAY RF Latex 100, • CORMAY RPR Carbon 500, • CORMAY TPHA 100. <p>The current labelling below:</p> <ul style="list-style-type: none"> • CORMAY CRP Latex 100: CRP-Latex contains boric acid. Danger H360FD May damage fertility. May damage the unborn child.  P201 Obtain special instructions before use. P308+P313 IF exposed or concerned: Get medical advice. P405 Store locked up. CRP-Control (+) and CRP-Control (-) contain 2-methyl-2H-isothiazol-3-one. Warning H317 May cause an allergic skin reaction. P280 Wear protective gloves, protective clothing, eye protection or face protection.  P302+P352 IF ON SKIN: Wash with plenty of soap and water. 		

	<ul style="list-style-type: none"> <p>• CORMAY ASO Latex 100: ASO-Control (+) and ASO-Control (-) contain 2-methyl-2H-isothiazol-3-one. Warning H317 May cause an allergic skin reaction. P280 Wear protective gloves, protective clothing, eye protection or face protection. P302+P352 IF ON SKIN: Wash with plenty of soap and water.</p>  <p>• CORMAY RF Latex 100: RF-Control (+) and RF-Control (-) contain 2-methyl-2H-isothiazol-3-one. Warning H317 May cause an allergic skin reaction. P280 Wear protective gloves, protective clothing, eye protection or face protection. P302+P352 IF ON SKIN: Wash with plenty of soap and water.</p>  <p>• CORMAY RPR Carbon 500: RPR-Control (+) and RPR-Control (-) contain 2-methyl-2H-isothiazol-3-one. Warning H317 May cause an allergic skin reaction. P280 Wear protective gloves, protective clothing, eye protection or face protection. P302+P352 IF ON SKIN: Wash with plenty of soap and water.</p>  <p>• CORMAY TPHA 100: TPHA-Control (+) and TPHA-Control (-) contain 2-methyl-2H-isothiazol-3-one. Warning H317 May cause an allergic skin reaction. P280 Wear protective gloves, protective clothing, eye protection or face protection. P302+P352 IF ON SKIN: Wash with plenty of soap and water.</p> 
Reason for change	<p>Fulfilling the requirements of the Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006 taking into account Adaptation to Technical Progress (ATP) to the CLP Regulation.</p>

Implementation	September 2023
Additional information	<p>Products with the new labelling will be available for sale beginning from September 2023. Products with the old labelling may also be available – while stocks last.</p> <p>Once you receive the notification, follow the new classification of the reagents. Use personal protective equipment when working.</p> <p>Updated material safety data sheets (MSDS) are available at: https://cormaydiagnostics.com/customer-support/documents</p>
Attachments	N/A

If there are any questions or concerns with respect to this change, please feel to contact your Area Sales Managers

Quality Assurance Coordinator