

**CHANGE NOTIFICATION**

Notification number	<b>PZ_2401_01</b>	Date	08.01.2024
Change number	043		
Summary	MYTHIC 18 VET software update.		
Affected products	Cat no	Product name	
	N/A	MYTHIC 18 VET	
Change category	<input type="checkbox"/> raw materials <input type="checkbox"/> quality control <input type="checkbox"/> packages <input type="checkbox"/> informations provided with the product <input type="checkbox"/> manufacturing <input checked="" type="checkbox"/> final product <input type="checkbox"/> other		
Change description	New autonomy calculation with added workload input parameter, added Autoclean and Autorinse counters, updated cycles to increase cleaner consumption reliability and new splash screen will be implemented in the new release of software.		
Reason for change	MYTHIC 18 VET software enhancement.		
Implementation	January 2024		
Additional information	The MYTHIC 18 VET has to be updated with software version 1.2.1-004 downloaded according to the attached instruction (Attachment No. 1).		
Attachments	Attachment No. 1 Instruction – how to download the software.		

If there are any questions or concerns with respect to this change, please feel to contact your International Service Engineers.

Quality Assurance Manager