

CHANGE NOTIFICATION

Notification number	PZ_2401_02	Date	15.01.2024
Change number	105		
Summary	MYTHIC 70 software update.		
Affected products	Cat no	Product name	
	N/A	MYTHIC 70	
Change category	<input type="checkbox"/> raw materials <input type="checkbox"/> quality control <input type="checkbox"/> packages <input type="checkbox"/> informations provided with the product <input type="checkbox"/> manufacturing <input checked="" type="checkbox"/> final product <input type="checkbox"/> other		
Change description	For the MYTHIC 70 analyzer, the function of hiding flags has been changed - the option to hide the flag "!" on result's printout has been implemented in the new release of software: V0.5.0-005. The flags "!" will still be visible on the result displayed on the analyzer.		
Reason for change	MYTHIC 70 software enhancement.		
Implementation	January 2024		
Additional information	The MYTHIC 70 has to be updated with software version V0.5.0-005 downloaded according to the attached instruction (Attachment No. 1).		
Attachments	Attachment No. 1 Instruction – how to download the software.		

If there are any questions or concerns with respect to this change, please fell to contact your Area Sales Manager.

Quality Assurance Manager